



Integrating the Healthcare Enterprise

IHE Integration Statement

**InDex® Archive V. 2.5
with AccessNET (StorComm)**

10 October 2003

Copyright 2003 by InSiteOne, Inc.

IHE Integration Statement			
Vendor	Product Name	Version	Date
InSiteOne, Inc.	InDex® Archive with AccessNET (StorCOMM)	V. 2.5	10 October 2003
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	<ul style="list-style-type: none"> Image Manager/Image Archive Performed Procedure Step Manager 	none	
Patient Information Reconciliation	<ul style="list-style-type: none"> Image Manager/Image Archive Performed Procedure Step Manager 	none	
Consistent Presentation	<ul style="list-style-type: none"> Image Manager/Image Archive Performed Procedure Step Manager 	none	
Access to Radiology Information	<ul style="list-style-type: none"> Image Manager/Image Archive Report Repository 	none	
Key Image Notes	<ul style="list-style-type: none"> Image Manager/Image Archive 	none	
Report Workflow	<ul style="list-style-type: none"> Image Manager/Image Archive 	none	
Simple Image and Numeric Report	<ul style="list-style-type: none"> Report Repository 	none	
Evidence Documents	<ul style="list-style-type: none"> Image Manager/Image Archive 	none	
IHE information at InSiteOne, Inc: www.insiteone.com/aboutindex_documents.htm			
IHE information at StorComm, Inc: www.storcomm.com/Products/IntegrationSolutions/IHE_IntegrationStatement.htm			
Links to Standards Conformance Statements for the Implementation			
DICOM www.insiteone.com/aboutindex_documents.htm			
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	