



Integrating the Healthcare Enterprise

*IHE Integration Statement*

**InDex® Archive V. 5.1.2**

3 March 2009

Copyright 2009 by InSite One, Inc.

## IHE Integration Statement

Vendor	Product Name	Version	Date
InSite One, Inc.	InDex® Archive	V. 5.1.2	3 March 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Access to Radiology Information	Image Manager/Image Archive Report Repository	none	
Cardiac Catheterization Workflow	Image Manager/Image Archive Performed Procedure Step Manager	<ul style="list-style-type: none"> <li>• Cardiac Cath</li> <li>• Intermittently Connected Modality</li> </ul>	
Consistent Presentation of Images	Image Manager/Image Archive	none	
Consistent Time	Time Client	none	
Echocardiography Workflow	Image Manager/Image Archive Performed Procedure Step Manager	<ul style="list-style-type: none"> <li>• Echocardiography</li> <li>• Intermittently Connected Modality</li> </ul>	
Evidence Documents	Image Manager/Image Archive	none	
Evidence Documents Cardiology	Image Manager/Image Archive	none	
Key Image Notes	Image Manager/Image Archive	none	
Mammography	Image Manager/Image Archive	none	
Nuclear Medicine Image	Image Manager/Image Archive	none	
Patient Information Reconciliation	Image Manager/Image Archive Performed Procedure Step Manager	none	
Report Workflow	Report Repository	none	
Scheduled Workflow	Image Manager/Image Archive Performed Procedure Step Manager	none	
Simple Image and Numeric Report	Report Repository	none	
<b>IHE information at InSiteOne, Inc:</b> <a href="http://www.insiteone.com/why-dicom.php">www.insiteone.com/why-dicom.php</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
DICOM <a href="http://www.insiteone.com/why-dicom.php">www.insiteone.com/why-dicom.php</a>			
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/e">www.jira-net.or.jp/e</a>	