

IHE Integration Statement

Vendor	Product Name	Version	Date
InSite One	InDex Archive	6.1	January 24, 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager/Archive	
	Report Repository	
Audit Trail and Node Authentication	Secure Node	Intermittently Connected Modality, Cardiac Cath
Cardiac Catheterization Workflow	Image Manager/Archive	
Cardiology Evidence Documents	Image Manager/Archive	
	Consistent Presentation of Images	
Consistent Time	Time Client	Intermittently Connected Modality, Echocardiography
Cross-Enterprise Clinical Documents Share (XDS.b)	Document Source	
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	
Echocardiography Workflow	Image Manager/Archive	
	Performed Procedure Step Manager	
Evidence Documents	Image Manager/Archive	
Key Image Note		
Mammography Image	Performed Procedure Step Manager	
NM Image		
Patient Information Reconciliation	Image Manager/Archive	
Scheduled Workflow		
Simple Image and Numeric Report	Performed Procedure Step Manager	
	Report Repository	

Internet address for vendor's IHE information: <http://insiteone.com/why-dicom.php>

Links to Standards Conformance Statements for the Implementation

HL7	http://insiteone.com/pdf/InDex5.1.3HL7Conformance.pdf
DICOM	http://insiteone.com/pdf/InDex5.1.3%20DICOMConformance.pdf

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--